

Immunofluorescence labeling

Solutions:

MTSB (1x)

50 mM PIPES 15,1 g (PIPES free acid, mw=302,4)
5 mM MgSO₄·7H₂O 1,23 g
5 mM EGTA 10 ml of 0,5 M stock solution
add solid KOH to dissolve PIPES first
final volume 1000 ml
pH=6,9 (use KOH to set)

EGTA (0,5 M stock solution)

3,8 g EGTA in 10 ml H₂O, add solid KOH until clear
make to 20 ml, pH=8

PBS

0,14 M NaCl 8,0 g
2,7 mM KCl 0,2 g
6,5 mM Na₂HPO₄·2H₂O 1,15 g
1,5 mM KH₂PO₄ 0,2 g
make with dH₂O up to 1000 ml, pH=7,3

Formaldehyde fixation

18 ml MTSB
2 ml 37% formaldehyde

Wax (PEG + 1-hexadecanol = 9+1 w/w)

melt PEG 400 distearate at 65°C, add 1-hexadecanol
stirr 3-4 hrs
pour wax into Aluminium containers, cool at RT

First antibody

Dilute in PBS, add 10 µl of BSA stock solution per 1 ml of PBS

Second antibody

Dilute in PBS, add 10 µl of BSA stock solution per 1 ml of PBS

BSA stock

100 mg in 1 ml dH₂O, store 100 µl aliquots at -20°C

DAPI stock (100 mM)

1 mg in 10 ml H₂O

DAPI working solution

2 µl of 100 mM stock solution in 2 ml PBS

Mounting medium

100 mg p-phenylenediamine in 10 ml PBS, add 90 ml of glycerol
set pH to 8,0 with carbonate/bicarbonate or Tris pH=8,0
store at -20°C

Procedure:

1. Formaldehyde fixation	60 min
2. MTSB	30 min
3. PBS	15 min
4. Dehydratation	
30 % ethanol in PBS	30 min
50 % ethanol in PBS	30 min
70 % ethanol in PBS	30 min
90 % ethanol in PBS	30 min
97 % ethanol	30 min
Toluidine blue in ethanol	10 min
5. Embedding	
97 % ethanol	10 min
wax + ethanol (1:1)	overnight
100 % wax	1-2 hrs
6. Sectioning	
sections attached to glycerol-albumine coated slides	overnight
7. Dewaxing and Rehydratation	
97 % ethanol	3x10 min
90 % ethanol in PBS	10 min
50 % ethanol in PBS	10 min
PBS	2x10 min
8. Labeling	
First antibody (at RT)	1 hr
PBS	
Second antibody (at RT)	1 hr
PBS	10 min
DAPI	10 min
PBS	10 min
Toluidine Blue (0,01 % in PBS)	10 min
PBS	10 min
9. Apply mounting medium and coverslide, seal with nail varnish	